INDEX, ÍNDICE, ΕΥΡΕΤΗΡΙΟ	
1-20	ENGLISH
21 -41	ESPAÑOL
42 -61	ITALIANO
62 -82	DEUTSCHE
83 - 102	CZECH
103 - 124	ΕΛΛΗΝΙΚΟ
125 - 145	FRANÇAIS
146 -166	POLSKI
167 - 185	中文

Instructions for Use

Thermo Scientific Sensititre™ YeastOne™ YO10 & ITAMYUCC Susceptibility Plates

<i>Sensititre™ YeastOne™ YO10 & ITAMYUCC</i>	
Revision Date:	<i>August 2024</i>

THERMO SCIENTIFIC SENSITITRE YEASTONE SUSCEPTIBILITY PLATES – YO10 & ITAMYUCC

For *in vitro* Diagnostic Use

Sensititre plates are 96-well microtitre plates, available in both standard and custom formats. Multiple antimicrobials/antifungals are available in various dilution ranges (on a wide variety of plate formats). Sensititre YeastOne Plates are used to determine MIC (Minimum Inhibitory Concentration).

After inoculation with a standardised suspension of organisms, the plate is covered with an adhesive seal, and incubated at 35°C for 24 hours. The MICs for the test organisms are read visually with respect to a colourimetric change in order to determine the lowest antifungal concentration that inhibits the growth of the isolate.

This instruction for use (IFU) document covers two plates from the Sensititre YeastOne range, YO10 and ITAMYUCC, which will be referred to as 'the plate'. Unless specifically stated, all information is applicable to both plates.

Intended Use

YO10 and ITAMYUCC are *in vitro* diagnostic devices intended for use with cultured samples from patients suspected of having yeast infection.

The plate is a semi-quantitative test for susceptibility testing of clinical isolates of *Candida* spp. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a yeast infection, the device is not automated and requires a pure, agar grown culture. This device is for professional use only and is not a companion diagnostic.

Principles of Method

The clinical specimen, once isolated and purified from the patient specimen, is inoculated at a known cell density into a growth medium and incubated for 24 hours in the presence of a suitable range of antifungal concentrations. The lowest concentration of the log dilution range that does not show colour change is considered to be the MIC. This MIC can then be linked to a positive therapeutic treatment, utilizing regulatory approved breakpoint concentrations, based on the identification of the organism and body site of the infection. The intended user can then use this information to decide the best course of action based on this and other patient specific information.

Materials Provided

- 10 x Sensititre plates
- 10 x Adhesive seals

Plate layouts

Please refer to figures 1 and 2 for antimicrobial positioning and concentration (µg/ml). For drug abbreviations please refer to table 1.

Figure 1: Drug layout and concentrations in Sensititre YO10

	1	2	3	4	5	6	7	8	9	10	11	12
A	POS	AND 0.015	AND 0.03	AND 0.06	AND 0.12	AND 0.25	AND 0.5	AND 1	AND 2	AND 4	AND 8	AB 0.12
B	MF 0.008	MF 0.015	MF 0.03	MF 0.06	MF 0.12	MF 0.25	MF 0.5	MF 1	MF 2	MF 4	MF 8	AB 0.25
C	CAS 0.008	CAS 0.015	CAS 0.03	CAS 0.06	CAS 0.12	CAS 0.25	CAS 0.5	CAS 1	CAS 2	CAS 4	CAS 8	AB 0.5
D	FC 0.06	FC 0.12	FC 0.25	FC 0.5	FC 1	FC 2	FC 4	FC 8	FC 16	FC 32	FC 64	AB 1
E	PZ 0.008	PZ 0.015	PZ 0.03	PZ 0.06	PZ 0.12	PZ 0.25	PZ 0.5	PZ 1	PZ 2	PZ 4	PZ 8	AB 2
F	VOR 0.008	VOR 0.015	VOR 0.03	VOR 0.06	VOR 0.12	VOR 0.25	VOR 0.5	VOR 1	VOR 2	VOR 4	VOR 8	AB 4
G	IZ 0.015	IZ 0.03	IZ 0.06	IZ 0.12	IZ 0.25	IZ 0.5	IZ 1	IZ 2	IZ 4	IZ 8	IZ 16	AB 8
H	FZ 0.12	FZ 0.25	FZ 0.5	FZ 1	FZ 2	FZ 4	FZ 8	FZ 16	FZ 32	FZ 64	FZ 128	FZ 256

Figure 2: Drug layout and concentrations in Sensititre ITAMYUCC

	1	2	3	4	5	6	7	8	9	10	11	12
A	POS	AND 0.015	AND 0.03	AND 0.06	AND 0.12	AND 0.25	AND 0.5	AND 1	AND 2	AND 4	AND 8	AB 0.12
B	MF 0.008	MF 0.015	MF 0.03	MF 0.06	MF 0.12	MF 0.25	MF 0.5	MF 1	MF 2	MF 4	MF 8	AB 0.25
C	CAS 0.008	CAS 0.015	CAS 0.03	CAS 0.06	CAS 0.12	CAS 0.25	CAS 0.5	CAS 1	CAS 2	CAS 4	CAS 8	AB 0.5
D	ISA 0.008	ISA 0.015	ISA 0.03	ISA 0.06	ISA 0.12	ISA 0.25	ISA 0.5	ISA 1	ISA 2	ISA 4	ISA 8	AB 1
E	PZ 0.008	PZ 0.015	PZ 0.03	PZ 0.06	PZ 0.12	PZ 0.25	PZ 0.5	PZ 1	PZ 2	PZ 4	PZ 8	AB 2
F	VOR 0.008	VOR 0.015	VOR 0.03	VOR 0.06	VOR 0.12	VOR 0.25	VOR 0.5	VOR 1	VOR 2	VOR 4	VOR 8	AB 4
G	IZ 0.015	IZ 0.03	IZ 0.06	IZ 0.12	IZ 0.25	IZ 0.5	IZ 1	IZ 2	IZ 4	IZ 8	IZ 16	AB 8
H	FZ 0.12	FZ 0.25	FZ 0.5	FZ 1	FZ 2	FZ 4	FZ 8	FZ 16	FZ 32	FZ 64	FZ 128	FZ 256

Table 1: Abbreviations of antimicrobics used in YO10 and ITAMYUCC

Abbreviation	Antimicrobics
POS	Positive Control
MF	Micafungin
CAS	Caspofungin
FC	5-Flucytosine
PZ	Posaconazole
VOR	Voriconazole
IZ	Itraconazole
FZ	Fluconazole
AND	Anidulafungin
AB	Amphotericin B
ISA	Isavuconazole

Materials Required But Not Provided

For a list of materials not provided with the device but required to correctly operate Sensititre YeastOne YO10 & ITAMYUCC plates, please refer to table 2.

Table2: *Description of materials not provided with Sensititre YeastOne plates*

Component	Product code	Description
Sensititre YeastOne Broth	Y3462	11ml Yeast broth for use with Sensititre YeastOne plates
Sensititre Demineralised Water	T3339	5ml demineralised water for use with Sensititre plates.
0.5 McFarland Turbidity Standard	E1041	Standard used for preparing bacterial suspensions
Bacteriological Loop	N/A	Inoculating loop for use in bacterial culture & preparation
20µl Pipette and Disposable Tips	N/A	A laboratory tool commonly used to transport and dispense a measured volume of liquid within a disposable tip
Sterile Inoculum Reservoir	N/A	Sterile reservoir to hold microbial suspensions
100µl Pipette and Disposable Tips	N/A	A laboratory tool commonly used to transport and dispense a measured volume of liquid within a disposable tip

Optional Instruments Not Provided

Instruments that can be used with Sensititre YeastOne plates, but are not essential, are listed below:

Table 3: *Optional Instruments that may be used with Sensititre YeastOne Plates*

Component	Product Code	Description
Sensititre Nephelometer	V3011	The Sensititre Nephelometer is used for standardizing the inoculum density of a fungal suspension. Light from an ultra-bright light emitting diode (LED) is scattered by particles or microbial cells in suspension and detected by a photodiode. The quantity of scattered light is compared with a reference value obtained by calibrating the Nephelometer with a 0.5 McFarland turbidity standard
Thermo Scientific Sensititre AIM™ (Automated Inoculation Delivery System)	V3020	The Sensititre AIM is a microprocessor-controlled instrument that rapidly delivers inoculum in multiples of 50µl to a 96-well microtitration plate.
Thermo Scientific Sensititre Dosing Heads	E3010	Disposable used with the AIM to deliver inoculum to each well of the Sensititre plate.
Manual Viewer	V4007	Mirrored manual viewbox for performing simple visual reads of Sensititre plates
Thermo Scientific Sensititre Digital MIC Viewing System	V2021	The Vizion is an instrument that takes a digital image of a Sensititre susceptibility plate and sends it to the SWIN software so that the results are displayed at a SWIN PC and can be read by the user. The Vizion is linked via a supplied cable to a USB port on a computer (PC) running the SWIN software.

Storage, Shelf Life and Handling Conditions

Plates should be stored at room temperature (15-25°C) away from direct sunlight and direct heat. Each plate is packaged in foil with a silica gel desiccant. Do not use the plate if (1) the desiccant colour is not orange, (2) the expiration date has passed, or (3) the foil pouch is damaged. Inoculate plate within 5 hours of removal from pouch.

Warning and Precautions

Results should be used as an aid in selecting the drug of choice for treatment. Only trained personnel should operate the system. Proper organism handling and disposal methods should be used. If the foil packaging is damaged or the desiccant is not orange, do not use the device.

When applying the adhesive seal, avoid creasing as this may cause erroneous result or “well skips”. A blue well in a series of pink growth wells indicates a skip/erroneous result and should be ignored. The MIC should be read above any skipped wells. If there is more than one skipped well, the antifungal should not be reported.

Follow incubation procedures listed below as incubation temperatures over 35°C may affect performance.

There are no reasonably foreseeable external influences that will impact the correct working of this device during routine use. The device is sensitive to humidity and temperature and should be used in accordance with the proper storage and handling conditions stated in this IFU and on the packaging.

There is some residual risk of user or environment being exposed to toxic/hazardous/ infectious materials. These risks are mitigated by trained clinical professionals using the device, following established precautions for microbial hazards and the warnings listed in the Safety Data Sheet (SDS) and IFU.

The device does not contain CMR, endocrine disrupting substances or materials likely to result in an allergic reaction under normal use.

There is no infectious material included with the device. Cultures, containers and other contaminated materials must be sterilised after use in accordance with guidelines for the handling and disposal of biohazardous waste.

Refer to the SDS available at www.trekds.com/techinfo for safe handling and disposal of the product.

Interfering substances and limitations

1. Sensititre YeastOne Plates YO10 and ITAMYUCC are for use with *Candida* species only. They are not intended for *Cryptococcus* species, *Aspergillus* species and miscellaneous rapid growing yeast species. They are not intended for fastidious or slow growing yeast such as *Histoplasma* species or *Blastomyces* species, and filamentous fungi.
2. Due to the low performance of fluconazole, anidulafungin and voriconazole with *C. tropicalis*, isolates of *C. tropicalis* should be tested with an alternate method.
Comparison between the Sensititre YeastOne System at 24 hours and the CLSI reference method at 48 hours was evaluated. However, due to the difficulty in correlating end points of trailing organisms (*C. albicans*) at 48 hours incubation, high error rates were observed. The CLSI official standards, M27-Ed4 2017 recommends that readings should be performed at 24 hours for clinical isolates and QC strains.
3. Testing of fungi and antifungal agents is inherently less precise than testing bacteria.
4. For additional guidance, refer to CLSI Antifungal Susceptibilities for yeast, Standard M27.
5. Colour change is the indicator of the end point, not turbidity. This feature helps address concerns with the interpretation of certain *Candida* species because of 'trailing'. Trailing is more commonly seen with isolates other than those of blood and other sterile body fluids.
6. Do not read at 24 hours if the control well has not undergone a complete colour change from blue to pink.
7. Use only with Sensititre System approved yeast susceptibility inoculum broth. The use of other broths could result in error.
8. As with any *in vitro* susceptibility testing method, the results of testing should be correlated with the patient's clinical response to prescribed therapy.
9. Correlation of the MIC for caspofungin to the treatment outcome following caspofungin use has not been fully established⁹.
10. Only instruments supported by the Sensititre system, i.e. a simple mirror viewer, Sensititre Vizion, must be used to report results with Sensititre products. Any other system used will not be supported.
11. Itraconazole can occasionally come out of solution at concentrations ≥ 4 µg/ml. This can result in the affected well exhibiting growth and turning pink.

Antifungal limitations

Antifungal Agent(s)	Organism	Limitation
Caspofungin	<i>Candida</i> spp. excluding <i>Candida glabrata</i>	The ability of the Sensititre plate containing caspofungin to determine resistant strains has not currently been established due to testing an insufficient number of resistant isolates.
Fluconazole	<i>Candida glabrata</i>	The ability of the Sensititre plate containing fluconazole to determine resistant strains has not currently been established due to testing an insufficient number of resistant isolates.
Caspofungin	<i>Candida tropicalis</i> and <i>Candida albicans</i>	MIC values with the Sensititre YeastOne Susceptibility System with caspofungin in the dilution range of 0.015-16 µg/mL tended to be in exact agreement or at least one doubling dilution higher for <i>C. tropicalis</i> , <i>C. albicans</i> , and <i>C. glabrata</i> compared to the CLSI broth microdilution reference method.
Anidulafungin Fluconazole Voriconazole	<i>Candida tropicalis</i>	Due to poor performance of the Sensititre YeastOne Susceptibility System with anidulafungin, fluconazole and voriconazole in the comparator study, <i>C. tropicalis</i> isolates should be tested with an alternative method.
Caspofungin	<i>Candida glabrata</i>	Due to categorical agreement below 90% with the Sensititre YeastOne Susceptibility System when compared to the CLSI antifungal broth microdilution method caused by the occurrence of categorical errors and to avoid potential false susceptible results, perform an alternative method of testing prior to reporting results for caspofungin and <i>C. glabrata</i> when the MIC is ≤ 0.12µg/ml.
Itraconazole	<i>Candida lusitanae</i>	Due to the low performance (EA <90%) of itraconazole with <i>C. lusitanae</i> , isolates of <i>C. lusitanae</i> should be tested with an alternative method.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Handling and preparation of the specimen

Specimens should be collected, transported, stored and plated onto primary isolation medium using standard procedures¹.

Inoculation Procedure

Allow all broths to come to room temperature before use. Plates should be inoculated within 5 hours of removal from the pouch. A final organism density of approximately $1.5 - 8 \times 10^3$ CFU/ml is recommended. Sensititre System approved broths are performance tested for use with the Sensititre System. The Sensitre YeastOne System is currently only suitable for MIC determinations using CLSI reference methods and standards.

Steps 1 and 2 should be completed within 15 minutes.

1. Pick several well isolated colonies of >1mm diameter from a pure 24-hour culture (Sabouraud Dextrose Agar) of the yeast isolate and emulsify into demineralised water. Mix well ensuring the suspension is uniform, vortex if required. If clumping occurs, allow the suspension to settle before adjusting the density. Adjust to a 0.5 McFarland standard visually or with a Sensititre Nephelometer.
2. Transfer 20 µl of the suspension into 11 ml of Sensititre YeastOne Broth to give a final inoculum of $1.5 - 8 \times 10^3$ CFU/ml.
3. Transfer 100µl of the final suspension to the Sensititre YeastOne Plate within 15 minutes of completing step 2 using either:
 - a. **Sensititre AIM-** Replace the tube cap with a single-use Sensititre dosing head and inoculate the plate according to the Sensititre AIM user manual. Remove the test tube dosing head combination from the Sensititre AIM System within 30 seconds of dosing a plate and store inverted in a rack or discard.
 - b. **Manual pipette –** Pour the broth into a sterile seed trough and inoculate the plate using an appropriate pipette.
4. A check of the colony count should be done by removing 10µl from the positive control well and plating onto Sabouraud Dextrose Agar (SDA). A correct inoculum will produce 15-80 colonies.
5. Cover all wells with the adhesive seal. Avoid creases when applying the adhesive seal.

Incubation

After inoculation, incubate plates at 35 °C in a non-CO₂ incubator for 24 to 25 hours. Incubation temperatures over 35°C may affect performance.

Interpretation of results and calculation of analytical result

Plates may be read visually under normal laboratory lighting, using a manual mirror viewer or by using the Sensititre Vizion System (refer to the Sensititre Vizion System user manual for additional information).

Yeast growth in the antifungal solutions will be evident as a change in the colorimetric growth indicator from blue (negative) to pink (positive). Some yeast species may not change the indicator completely to pink but display an intermediate purple colour. Some organisms may exhibit purple colours with fluconazole, itraconazole, ketoconazole, posaconazole and voriconazole. A colour change to a bleached pink or purple is associated with good growth and should be read as positive.

To interpret results:

1. Examine the positive growth well after incubation.
2. If the positive growth well is pink, the endpoints for the antifungals can be determined. If the well is blue or faint purple, re-incubate for an additional 24 hours and re-examine.

Do not read turbidity in the Sensititre YeastOne plate, read only colour change.

3. The MIC is the lowest concentration of an antifungal agent that substantially inhibits growth of the organism as detected by a colour change from pink to blue.

Interpretation of Results

Table 4. Examples of typical test results and how to interpret them.

	Well Concentration µg/ml						P = PINK: Positive growth indication B = BLUE: Negative growth indication
	1	2	4	8	16	32	
A.	P	P	P	B	B	B	Typical growth pattern; MIC endpoint is 8 µg/ml.
B.	P	P	P	P	P	P	Growth in all wells; MIC endpoint is >32 µg/ml.
C.	B	B	B	B	B	B	No growth in any well; MIC endpoint is ≤1 µg/ml.
D.	P	P	P	B	P	P	"Skipped Well". MIC endpoint is >32 µg/ml. Disregard "skip" when wells on either side have growth. If more than one "skip" should occur in a column, the test results are invalidated*
E.	P	P	B	B	P	P	Double "Skipped Well". The test should be repeated *

*With careful technique, these occurrences are not common.

If MICs for *Candida* spp. are measured using a scale that yields results falling between categories, the next higher category is implied. Thus, an isolate with a fluconazole MIC of 12.5 µg/mL would be placed in the S-DD category (19).

Reading Notes

Amphotericin B.

For amphotericin B at 24 hours, the endpoints are typically well defined, and the MIC is read as the lowest drug concentration that prevents any discernible colour change. Trailing endpoints with Amphotericin B are not usually encountered.

The first well showing a distinct colour change as compared to the positive growth well is the MIC.

Flucytosine and Azole Antifungals

Candida albicans, *C. glabrata* and *C. tropicalis* with flucytosine and azoles, such as fluconazole, itraconazole, ketoconazole, voriconazole and posaconazole may give endpoints that are typically less sharp because of trailing growth, and may be a significant source of variability. Trailing occurs when a slight colour change persists and it is often identical for all drug concentrations above the MIC. The MIC should be read as the first well showing a less intense colour change (indicating a 50% decrease in growth) compared to the positive growth control well.

Reference strains of defined susceptibility may also help to train personnel. Isolates of *Candida krusei* are assumed to be intrinsically resistant to fluconazole and their MICs should not be interpreted¹. A comment should accompany the test result reported.

Trailing endpoint: This occurs when a slight colour change persists and is often identical in several concentrations. The MIC should be read as the first well showing a less intense colour change compared to the positive growth control well (well after a complete pink colour change).

Echinocandins

MIC end points should be determined after 24 hours of incubation at 35°C. The MIC should be read as the first well showing a less intense colour change as compared to the positive growth control well.


Itraconazole

Itraconazole can occasionally precipitate at concentrations ≥ 4 $\mu\text{g/ml}$ (well highlighted with the black square). This can result in the affected well exhibiting growth and turning pink.


Occasionally unexpected growth in the higher concentrations of Itraconazole on the Sensititre YeastOne plates may be observed which results in pinking of these wells.

This paradoxical effect, known as the Eagle phenomenon¹⁸, occurs when an increase in the antimicrobial concentration results in an increase in the number of micro-organisms that survive. This may be due a self-antagonising effect with the associated binding receptor (e.g. penicillin binding proteins in the case of a penicillin).

Resolution

The growth in the high concentration should be ignored unless you have growth in all of the other concentrations of Itraconazole. In the example below the well highlighted with the black square is where the MIC result should be recorded.

If you have any other questions or concerns, please contact the technical support department on Phone: +44 (0) 1256 694287 | Fax: +44 (0) 1256 463388 or below for a list of contact email addresses.

Contamination & Skips

Note that a pink (growth) well between blue (no growth) wells may be indicative of contamination. Sub-culture well contents to ascertain the cause.

A blue well in a series of pink growth wells indicates a “skip” and should be ignored. The MIC should be read above any skip wells. If there is more than one skipped well, the antifungal should not be reported.

Quality Control Procedures

Frequency of quality control testing should be conducted according to local guidelines¹.

Inoculum should be cultured onto a suitable medium to check for purity. Test results are invalid if a mixed culture is detected.

All Sensititre plates include positive control wells. Tests are invalid unless there is distinct growth in the positive control well.

Control Materials/Strains

The following cultures from the American Type Culture Collection (ATCC™) are recommended for user quality control:

- *Issatchenkia orientalis* (*Candida krusei*) - ATCC 6258
- *Candida parapsilosis* - ATCC 22019

Clinical results should **not** be reported if QC results are not in the expected range.

The inoculation, reading and interpretation of Sensititre YeastOne susceptibility plates when testing for user quality control should be performed as described in the preceding section.

Table 5. Recommended 24 and 48 hour MIC limits for two quality control strains as per Broth Microdilution CLSI M27². Ranges that are different from or additional to published quality control ranges are underlined.

Antifungal Agent	<i>Issatchenkia orientalis</i> ATCC6258™		<i>Candida parapsilosis</i> ATCC 22019™	
	24 hours	48 hours	24 hours	48 hours
5-Flucytosine	4-16	8-32	<u>0.12-0.5*</u>	0.12-0.5
Amphotericin B	0.5-2	1-4	0.25-2	0.5-4
Anidulafungin	0.03-0.12	-	0.25-2	-
Caspofungin	0.12-1	0.25-1	0.25-1	0.5-4
Fluconazole	8-64	16-128	0.5-4	2-8
Itraconazole	0.12-1	0.25-1	0.06-0.5	0.06-0.5
Micafungin	0.06-0.5*	0.12 – 0.5	0.5-2	0.5-4
Posaconazole	0.06-0.5	0.12-1	0.03-0.25	0.06-0.25
Voriconazole	0.06-0.5	0.12-1	0.015-0.12	0.03-0.25
Isavuconazole	0.06-0.5	0.12-1.0*	0.015-0.06	0.03-0.12*

* Sensititre range (an internally generated QC range generated at the point of product development using CLSI M27 methodology)

Expected QC values are provided in CLSI M27M44S¹⁹.

Please refer to CLSI¹ for more information concerning interpretation of results.

Specific Performance Characteristics

Panels read manually are designed to give comparable performance to CLSI reference microbroth procedure. Comparable performance is defined as > 90% essential and categorical agreement to within two doubling dilutions of the reference MIC¹.

Precision of Measurement:

Essential agreement when testing QC organisms on Sensititre plates can be found in table 6. Essential agreement is reported when results are within 2 wells (serial dilution) of the model MIC or frozen reference result.

Table 6: *Essential agreement of QC organisms in Sensititre plates*

Antimicrobial	Number of tests	Results showing Essential Agreement
Amphotericin B	126	126 (100%)
Anidulafungin	154	154 (100%)
Caspofungin	32	32 (100%)
Fluconazole	123	123 (100%)
5-Flucytosine	155	155 (100%)
Itraconazole	135	135 (100%)
Micafungin	22	22 (100%)
Posaconazole	32	32 (100%)
Voriconazole	40	40 (100%)
Isavuconazole	144	144(100%)

Repeatability & Reproducibility:

Repeat daily testing conducted over three separate time periods is summarised in table 7. Repeatability and reproducibility of Sensititre plates is ≥95% when comparing to the CLSI M27M44S QC ranges or in-house Sensititre QC ranges.

Table 7: Manual data from Sensititre plates over three separate time periods

Antimicrobial	Number of repeats	Results falling within the CLSI M27M44S or in-house Sensititre QC ranges.
Amphotericin B	67	67 (100%)
Anidulafungin	35	35 (100%)
Caspofungin	112	112 (100%)
Fluconazole	53	53 (100%)
5-Flucytosine	35	35 (100%)
Itraconazole	35	35 (100%)
Micafungin	30	30 (100%)
Posaconazole	96	96 (100%)
Voriconazole	35	35 (100%)
Isavuconazole	144	144 (100%)

Analytical sensitivity

Manual susceptibility tests performed with QC organisms on Sensititre plates show ≥95% agreement with CLSI M27M44S QC ranges or in-house Sensititre ranges.

Table 8: Results from manual susceptibility testing performed on QC organisms.

Antimicrobial	Number of Tests	Results falling within the CLSI M27M44S or in-house Sensititre QC ranges.
Amphotericin B	126	126 (100%)
Anidulafungin	154	154 (100%)
Caspofungin	112	112 (100%)
Fluconazole	123	123 (100%)
5-Flucytosine	155	155 (100%)
Itraconazole	135	135 (100%)
Micafungin	100	100 (100%)
Posaconazole	112	112 (100%)
Voriconazole	115	115 (100%)
Isavuconazole	144	144 (100%)

Analytical Specificity

Analytical specificity is not applicable for this device as a pure culture is used and there is no cross reactivity and/or endogenous/exogenous interference sources.

Limit of Detection

The limit of detection for the plate is 1 log dilution in µg/ml (e.g. 0.03 – 0.06). Anything in between these concentrations would not be detected.

Table 9: Clinical performance characteristics for the Sensititre plates containing the following antifungals assessed using the CLSI reference method and breakpoints.

Antimicrobial Agent	Number of Isolates Tested	Category Agreement CLSI	Category Agreement (%) CLSI	VME %	ME %	mE%
Amphotericin B	N/A*	N/A*	N/A*	N/A*	N/A*	N/A*
Anidulafungin ^a	449	424	99.3%	0.0	0.0	3.0
Caspofungin	573	572	99.8%	0.0	1.0	0.0
Fluconazole ^a	1019	667	92.6%	1.89	0.92	3.53
5-Flucytosine	N/A*	N/A*	N/A*	N/A*	N/A*	N/A*
Itraconazole	N/A*	N/A*	N/A*	N/A*	N/A*	N/A*
Micafungin	315	314	99.7	0.0	0.0	0.3
Posaconazole	N/A*	N/A*	N/A*	N/A*	N/A*	N/A*
Voriconazole ^a	298	204	94.7%	5.0 ^b	0.0	4.7
Isavuconazole	N/A*	N/A*	N/A*	N/A*	N/A*	N/A*

*CLSI M60 (now M27M44S) 2nd edition breakpoints (BP) are no longer available.

a – Due to the low performance of fluconazole, voriconazole and anidulafungin with *C. tropicalis*, isolates of *C. tropicalis* should be tested with an alternate method.

b – There was one very major error for *C. albicans* which was considered to be related to incorrect interpretation of the reference method because of trailing endpoints and was therefore considered a random error.

c- The performance characteristics of the stated antifungals was only assessed against CLSI regulatory criteria. The data available was insufficient or incompatible to determine compliance with the EUCAST method.

Another source of clinical data was also collected from 2015-2021 with the performance characteristics described in table 10. below. This was initially conducted to assess the antifungals on Sensititre plates against EUCAST methodology, however changes to EUCAST guidelines and methodologies regarding tissue-treated plates resulted in the data being incompatible for this purpose. The data obtained were also assessed using the CLSI methodology and reference method. This assessment is presented below.

Table 10. Clinical performance characteristics for the Sensititre plates containing the following antifungals obtained from a 2015-2021 assessment using the CLSI reference method and breakpoints.

Antimicrobial Agent	Number of Isolates Tested	Category Agreement	Category Agreement (%)	Essential Agreement	Essential Agreement (%)
Anidulafungin ^a	1150	940	95.6	1072	93.2
Caspofungin	499	401	93.9	493	98.8
Fluconazole ^a	1150	727	91.9	1116	97.0
5-Flucytosine	498	N/A*	N/A*	493	99.0
Itraconazole	1148	N/A*	N/A*	1082	94.3
Micafungin	651	545	98.0	646	99.2
Posaconazole	499	N/A*	N/A*	481	96.4
Voriconazole ^a	1151	674	89.2 ^b	1112	96.6
Isavuconazole	499	N/A*	N/A*	484	97.0

*CLSI M60 (now M27M44S) 2nd edition breakpoints (BP) are no longer available.

a – Due to the low performance of fluconazole, voriconazole and anidulafungin with *C. tropicalis*, isolates of *C. tropicalis* should be tested with an alternate method.

b – The category agreement for voriconazole falls just under the >90% CLSI criterion due to *C. tropicalis* isolates performing poorly (74.3%), reducing the overall CA for this antifungal. As *C. tropicalis* isolates have a limitation for testing with voriconazole this CA can be justified.

Manufacturer Contact Information

The latest information can be obtained from www.trekds.com/techinfo or by contacting Thermo Fisher Scientific Microbiology Technical Services.

TREK Diagnostic Systems Ltd
Units 17-19,
Birches Industrial Estate,
East Grinstead,
West Sussex RH19 1XZ,
United Kingdom
Tel: +44 13423 18777

The ATCC Licensed Derivative® Emblem, the ATCC Licensed Derivative® word mark, and the ATCC catalog marks are trademarks of ATCC. Thermo Fisher Scientific Inc. is licensed to use these trademarks and to sell products derived from ATCC® cultures.

© 2021 Thermo Fisher Scientific Inc. All rights reserved. ATCC is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.


2797

Technical Support email addresses

Austria	microbiology.techsupport.de@thermofisher.com
Belgium	microbiology.techsupport.benelux@thermofisher.com
Czech Republic	mikrobiologie.tech.podpora.cz@thermofisher.com
Denmark	microbiology.techsupport.nordic@thermofisher.com
France	microbiologie.techsupport.fr@thermofisher.com
Finland	mikrobiologia.tekninentuki.fi@thermofisher.com
Germany	microbiology.techsupport.de@thermofisher.com
Italy	microbiologia.supportotecnico.it@thermofisher.com
Netherlands	microbiology.techsupport.benelux@thermofisher.com
Norway	microbiology.techsupport.nordic@thermofisher.com
Spain	microbiologia.soporte.es@thermofisher.com
Sweden	microbiology.techsupport.nordic@thermofisher.com
UK Technical Support	microbiology.techsupport.uk@thermofisher.com

Glossary of symbols

Symbol/Label	Meaning
	Manufacturer
	<i>In Vitro</i> Diagnostic Medical Device
	Temperature limitation (storage temp.)
	Lot / Batch Number
	Catalogue Number
	Do not re-use
	Consult Instructions for Use
	Contains sufficient for 10 tests
	Expiry Date
	Period after opening
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	Do not use if package is damaged. Consult instructions for use.
	Keep away from sunlight
	Authorized Representative in the European Community
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

Bibliography

- (1) CLSI M27. Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts (2017 or latest edition). CLSI (Clinical & Laboratory Standards Institute), 940 West Valley Road, Suite 1400, Wayne, PA 19087, US.
- (2) Espinel-Ingroff A *et al.* (1999). Multicentre comparison of the Sensititre YeastOne® Colorimetric Antifungal Plate with the NCCLS, M27-A Reference Method for Testing Clinical Isolates of Common and Emerging *Candida* spp., and other Yeasts and Yeast-like Organisms. *Journal of Clinical Microbiology*. **37**: 591-595.
- (3) Espinel-Ingroff A *et al.* (2002). Sensititre YeastOne® colorimetric antifungal panel for testing voriconazole against Isolates of *Candida* spp., a comparison with the NCCLS M27-A microdilution reference method. *Clinical Microbiology and Infection*. **8**: Suppl 1. Abstract 0125.
- (4) Davey K G *et al.* (1998). Comparison of a new commercial colorimetric microdilution method with a standard method for in-vitro susceptibility testing of *Candida* spp. and *Cryptococcus neoformans*., *Journal of Antimicrobial Chemotherapy*. **42**: 439-444.
- (5) Linares M *et al.* (2005). Susceptibility of filamentous fungi to voriconazole tested by two microdilution methods. *Journal of Clinical Microbiology*. **43**: 250-253.
- (6) Kartsonis N *et al.* (2005). Caspofungin susceptibility testing of isolates from patients with esophageal candidiasis or invasive candidiasis: relationship of MIC to treatment outcome. *Antimicrobial Agents and Chemotherapy*. **49**: 3616-3623.
- (7) Pfaller M A *et al.* (2004). Clinical evaluation of a dried commercially prepared microdilution panel for antifungal susceptibility testing of five antifungal agents against *Candida* spp. and *Cryptococcus neoformans*. *Diagnostic Microbiology and Infectious Diseases*. **50** :113-117.
- (8) Pfaller M A *et al.* (2004). Clinical evaluation of the Sensititre YeastOne® antifungal susceptibility testing of the new triazoles voriconazole, posiconazole and ravuconazole. *Journal of Clinical Microbiology*. **42**: 4577-4580.
- (9) Espinel-Ingroff A *et al.* (2004). Multicenter comparison of the Sensititre YeastOne® colorimetric antifungal panel with the NCCLS, M27-A2 reference method for testing new antifungal agents against clinical isolates of *Candida* spp. *Journal of Clinical Microbiology* **42**: 718- 721.
- (10) Linares M *et al.* (2004). Comparison of two microdilution methods for testing susceptibility of *Candida* spp. to voriconazole. *Journal of Clinical Microbiology* **42**: 899-902.
- (11) Espinel-Ingroff A *et al.* (1999). Multicenter comparison of the Sensititre YeastOne® colorimetric antifungal panel with the National Committee for Clinical Laboratory Standards M27-A reference method for testing clinical isolates of common and emerging *Candida* spp., *Cryptococcus* spp., and other Yeasts and Yeast-Like Organisms. *Journal of Clinical Microbiology* **37**: 591-595.
- (12) Barry A L *et al.* (2000). Quality control limits for broth microdilution susceptibility tests of ten antifungal agents. *Journal of Clinical Microbiology* **38**: 3457-3459.
- (13) Canton E *et al.* (2005). Sensititre YeastOne® Caspofungin susceptibility testing of *Candida* clinical isolates: Correlation with results of NCCLS, M27-A2 multicenter study . *Antimicrobial Agents and Chemotherapy*. **49**: 1604 - 1607.
- (14) Holliday N M *et al.* (2007) A reproducibility study with two echinocandins using the Sensititre YeastOne® susceptibility plate. *American Microbiology Abstract* **C-191**.
- (15) Torres-Nabona M *et al.* (2007) *In Vitro* activities of amphotericin B, caspofungin, itraconazole, posaconazole, and voriconazole against 45 clinical isolates of Zygomycetes: Comparison of NCCLS, M38-A, Sensititre_YeastOne® and the Etest. *Antimicrobial Agents and Chemotherapy*. **51**: 1126 -1129.

- (16) Espinel-Infgróff A (2006). Comparison of three commercial assays and a modified disk diffusion assay with two broth microdilution reference assays for testing Zygomycetes, *Aspergillus* spp, *Candida* spp, and *Cryptococcus neoformans* with posaconazole and Amphotericin B. *Journal of Clinical Microbiology* **44**: 3616-3622
- (17) Alexander B D *et al.* (2007). Comparative evaluation of Etest and Sensititre YeastOne® panels against the Clinical Laboratory Standards Institute M27-A2 reference broth microdilution method for testing *Candida* susceptibility to seven antifungal agents. *Journal of Clinical Microbiology* **45**: 698-706.
- (18) Eagle H. (1948). A paradoxical zone phenomenon in the bactericidal action of penicillin *in vitro*. *Science*, 1948; 107: 44-45
- (19) CLSI M27M44S. Performance Standards For Antifungal Susceptibility Testing Of Yeasts – 3rd Edition. CLSI (Clinical & Laboratory Standards Institute), 940 West Valley Road, Suite 1400, Wayne, PA 19087, US.

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***Placas de sensibilidad Sensititre™
YeastOne™YO10 e ITAMYUCC de Thermo Scientific***

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Made in the
United Kingdom

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- (1) CLSI M27. Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts (2017 or latest edition). CLSI (Clinical & Laboratory Standards Institute), 940 West Valley Road, Suite 1400, Wayne, PA 19087, US.
- (2) Espinel-Ingroff A *et al.* (1999). Multicentre comparison of the Sensititre YeastOne® Colorimetric Antifungal Plate with the NCCLS, M27-A Reference Method for Testing Clinical Isolates of Common and Emerging *Candida* spp., and other Yeasts and Yeast-like Organisms. *Journal of Clinical Microbiology*. **37**: 591-595.
- (3) Espinel-Ingroff A *et al.* (2002). Sensititre YeastOne® colorimetric antifungal panel for testing voriconazole against Isolates of *Candida* spp., a comparison with the NCCLS M27-A microdilution reference method. *Clinical Microbiology and Infection*. **8**: Suppl 1. Abstract 0125.
- (4) Davey K G *et al.* (1998). Comparison of a new commercial colorimetric microdilution method with a standard method for in-vitro susceptibility testing of *Candida* spp. and *Cryptococcus neoformans*., *Journal of Antimicrobial Chemotherapy*. **42**: 439-444.
- (5) Linares M *et al.* (2005). Susceptibility of filamentous fungi to voriconazole tested by two microdilution methods. *Journal of Clinical Microbiology*. **43**: 250-253.
- (6) Kartsonis N *et al.* (2005). Caspofungin susceptibility testing of isolates from patients with esophageal candidiasis or invasive candidiasis: relationship of MIC to treatment outcome. *Antimicrobial Agents and Chemotherapy*. **49**: 3616-3623.
- (7) Pfaller M A *et al.* (2004). Clinical evaluation of a dried commercially prepared microdilution panel for antifungal susceptibility testing of five antifungal agents against *Candida* spp. and *Cryptococcus neoformans*. *Diagnostic Microbiology and Infectious Diseases*. **50** :113-117.
- (8) Pfaller M A *et al.* (2004). Clinical evaluation of the Sensititre YeastOne® antifungal susceptibility testing of the new triazoles voriconazole, posiconazole and ravuconazole. *Journal of Clinical Microbiology*. **42**: 4577-4580.
- (9) Espinel-Ingroff A *et al.* (2004). Multicenter comparison of the Sensititre YeastOne® colorimetric antifungal panel with the NCCLS, M27-A2 reference method for testing new antifungal agents against clinical isolates of *Candida* spp. *Journal of Clinical Microbiology* **42**: 718- 721.
- (10) Linares M *et al.* (2004). Comparison of two microdilution methods for testing susceptibility of *Candida* spp. to voriconazole. *Journal of Clinical Microbiology* **42**: 899-902.
- (11) Espinel-Ingroff A *et al.* (1999). Multicenter comparison of the Sensititre YeastOne® colorimetric antifungal panel with the National Committee for Clinical Laboratory Standards M27-A reference method for testing clinical isolates of common and emerging *Candida* spp., *Cryptococcus* spp., and other Yeasts and Yeast-Like Organisms. *Journal of Clinical Microbiology* **37**: 591-595.
- (12) Barry A L *et al.* (2000). Quality control limits for broth microdilution susceptibility tests of ten antifungal agents. *Journal of Clinical Microbiology* **38**: 3457-3459.
- (13) Canton E *et al.* (2005). Sensititre YeastOne® Caspofungin susceptibility testing of *Candida* clinical isolates: Correlation with results of NCCLS, M27-A2 multicenter study. *Antimicrobial Agents and Chemotherapy*. **49**: 1604 - 1607.
- (14) Holliday N M *et al.* (2007) A reproducibility study with two echinocandins using the Sensititre YeastOne® susceptibility plate. *American Microbiology Abstract* **C-191**.

- (15) Torres-Nabona M *et al.* (2007) *In Vitro* activities of amphotericin B, caspofungin, itraconazole, posaconazole, and voriconazole against 45 clinical isolates of Zygomycetes: Comparison of NCCLS, M38-A, Sensititre_YeastOne® and the Etest. *Antimicrobial Agents and Chemotherapy*. **51**: 1126 -1129.
- (16) Espinel-Infgróff A (2006). Comparison of three commercial assays and a modified disk diffusion assay with two broth microdilution reference assays for testing Zygomycetes, *Aspergillus* spp, *Candida* spp, and *Cryptococcus neoformans* with posaconazole and Amphotericin B. *Journal of Clinical Microbiology* **44**: 3616-3622
- (17) Alexander B D *et al.* (2007). Comparative evaluation of Etest and Sensititre YeastOne® panels against the Clinical Laboratory Standards Institute M27-A2 reference broth microdilution method for testing *Candida* susceptibility to seven antifungal agents. *Journal of Clinical Microbiology* **45**: 698-706.
- (18) Eagle H. (1948). A paradoxical zone phenomenon in the bactericidal action of penicillin *in vitro*. *Science*, 1948; 107: 44-45
- (19) CLSI M27M44S. Performance Standards For Antifungal Susceptibility Testing Of Yeasts – 3rd Edition. CLSI (Clinical & Laboratory Standards Institute), 940 West Valley Road, Suite 1400, Wayne, PA 19087, US.

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TREK Diagnostic Systems Ltd
Units 17-19,
Birches Industrial Estate,
East Grinstead,
West Sussex RH19 1XZ,
Regno Unito
Tel: +44 (13423) 18777

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
Bibliografia

- (1) CLSI M27. Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts (2017 or latest edition). CLSI (Clinical & Laboratory Standards Institute), 940 West Valley Road, Suite 1400, Wayne, PA 19087, US.
- (2) Espinel-Ingroff A *et al.* (1999). Multicentre comparison of the Sensititre YeastOne® Colorimetric Antifungal Plate with the NCCLS, M27-A Reference Method for Testing Clinical Isolates of Common and Emerging *Candida* spp., and other Yeasts and Yeast-like Organisms. *Journal of Clinical Microbiology*. **37**: 591-595.
- (3) Espinel-Ingroff A *et al.* (2002). Sensititre YeastOne® colorimetric antifungal panel for testing voriconazole against Isolates of *Candida* spp., a comparison with the NCCLS M27-A microdilution reference method. *Clinical Microbiology and Infection*. **8**: Suppl 1. Abstract 0125.
- (4) Davey K G *et al.* (1998). Comparison of a new commercial colorimetric microdilution method with a standard method for in-vitro susceptibility testing of *Candida* spp. and *Cryptococcus neoformans*., *Journal of Antimicrobial Chemotherapy*. **42**: 439-444.
- (5) Linares M *et al.* (2005). Susceptibility of filamentous fungi to voriconazole tested by two microdilution methods. *Journal of Clinical Microbiology*. **43**: 250-253.
- (6) Kartsonis N *et al.* (2005). Caspofungin susceptibility testing of isolates from patients with esophageal candidiasis or invasive candidiasis: relationship of MIC to treatment outcome. *Antimicrobial Agents and Chemotherapy*. **49**: 3616-3623.
- (7) Pfaller M A *et al.* (2004). Clinical evaluation of a dried commercially prepared microdilution panel for antifungal susceptibility testing of five antifungal agents against *Candida* spp. and *Cryptococcus neoformans*. *Diagnostic Microbiology and Infectious Diseases*. **50**: 113-117.
- (8) Pfaller M A *et al.* (2004). Clinical evaluation of the Sensititre YeastOne® antifungal susceptibility testing of the new triazoles voriconazole, posiconazole and ravuconazole. *Journal of Clinical Microbiology*. **42**: 4577-4580.
- (9) Espinel-Ingroff A *et al.* (2004). Multicenter comparison of the Sensititre YeastOne® colorimetric antifungal panel with the NCCLS, M27-A2 reference method for testing new antifungal agents against clinical isolates of *Candida* spp. *Journal of Clinical Microbiology*. **42**: 718- 721.
- (10) Linares M *et al.* (2004). Comparison of two microdilution methods for testing susceptibility of *Candida* spp. to voriconazole. *Journal of Clinical Microbiology*. **42**: 899-902.
- (11) Espinel-Ingroff A *et al.* (1999). Multicenter comparison of the Sensititre YeastOne® colorimetric antifungal panel with the National Committee for Clinical Laboratory Standards M27-A reference method for testing clinical isolates of common and emerging *Candida* spp., *Cryptococcus* spp., and other Yeasts and Yeast-Like Organisms. *Journal of Clinical Microbiology*. **37**: 591-595.
- (12) Barry A L *et al.* (2000). Quality control limits for broth microdilution susceptibility tests of ten antifungal agents. *Journal of Clinical Microbiology*. **38**: 3457-3459.
- (13) Canton E *et al.* (2005). Sensititre YeastOne® Caspofungin susceptibility testing of *Candida* clinical isolates: Correlation with results of NCCLS, M27-A2 multicenter study . *Antimicrobial Agents and Chemotherapy*. **49**: 1604 -1607.
- (14) Holliday N M *et al.* (2007) A reproducibility study with two echinocandins using the Sensititre YeastOne® susceptibility plate. *American Microbiology Abstract* **C-191**.
- (15) Torres-Nabona M *et al.* (2007) *In Vitro* activities of amphotericin B, caspofungin, itraconazole, posaconazole, and voriconazole against 45 clinical isolates of Zygomycetes: Comparison of NCCLS, M38-A, Sensititre_YeastOne® and the Etest. *Antimicrobial Agents and Chemotherapy*. **51**: 1126 -1129.
- (16) Espinel-Ingroff A (2006). Comparison of three commercial assays and a modified disk diffusion assay with two broth microdilution reference assays for testing Zygomycetes, *Aspergillus* spp, *Candida* spp, and *Cryptococcus neoformans* with posaconazole and Amphotericin B. *Journal of Clinical Microbiology*. **44**: 3616-3622
- (17) Alexander B D *et al.* (2007). Comparative evaluation of Etest and Sensititre YeastOne® panels against the Clinical Laboratory Standards Institute M27-A2 reference broth microdilution method for testing *Candida* susceptibility to seven antifungal agents. *Journal of Clinical Microbiology*. **45**: 698-706.

- (18) Eagle H. (1948). A paradoxical zone phenomenon in the bactericidal action of penicillin *in vitro*. Science, 1948; 107: 44-45
- (19) CLSI M27M44S. Performance Standards For Antifungal Susceptibility Testing Of Yeasts – 3rd Edition. CLSI (Clinical & Laboratory Standards Institute), 940 West Valley Road, Suite 1400, Wayne, PA 19087, US.

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(15) Torres-Nabona M *et al.* (2007) *In Vitro* activities of amphotericin B, caspofungin, itraconazole, posaconazole, and voriconazole against 45 clinical isolates of Zygomycetes: Comparison of NCCLS, M38-A, Sensititre_YeastOne® and the Etest. *Antimicrobial Agents and Chemotherapy*. **51**: 1126 -1129.

(16) Espinel-Infgróff A (2006). Comparison of three commercial assays and a modified disk diffusion assay with two broth microdilution reference assays for testing Zygomycetes, *Aspergillus* spp, *Candida* spp, and *Cryptococcus neoformans* with posaconazole and Amphotericin B. *Journal of Clinical Microbiology* **44**: 3616-3622

(17) Alexander B D *et al.* (2007). Comparative evaluation of Etest and Sensititre YeastOne® panels against the Clinical Laboratory Standards Institute M27-A2 reference broth microdilution method for testing *Candida* susceptibility to seven antifungal agents. *Journal of Clinical Microbiology* **45**: 698-706.

(18) Eagle H. (1948). A paradoxical zone phenomenon in the bactericidal action of penicillin *in vitro*. *Science*, 1948; 107: 44-45

(19) CLSI M27M44S. Performance Standards For Antifungal Susceptibility Testing Of Yeasts – 3rd Edition. CLSI (Clinical & Laboratory Standards Institute), 940 West Valley Road, Suite 1400, Wayne, PA 19087, US.

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- (1) CLSI M27. Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts (2017 or latest edition). CLSI (Clinical & Laboratory Standards Institute), 940 West Valley Road, Suite 1400, Wayne, PA 19087, US.
- (2) Espinel-Ingroff A *et al.* (1999). Multicentre comparison of the Sensititre YeastOne® Colorimetric Antifungal Plate with the NCCLS, M27-A Reference Method for Testing Clinical Isolates of Common and Emerging *Candida* spp., and other Yeasts and Yeast-like Organisms. *Journal of Clinical Microbiology*. **37**: 591-595.
- (3) Espinel-Ingroff A *et al.* (2002). Sensititre YeastOne® colorimetric antifungal panel for testing voriconazole against Isolates of *Candida* spp., a comparison with the NCCLS M27-A microdilution reference method. *Clinical Microbiology and Infection*. **8**: Suppl 1. Abstract 0125.
- (4) Davey K G *et al.* (1998). Comparison of a new commercial colorimetric microdilution method with a standard method for in-vitro susceptibility testing of *Candida* spp. and *Cryptococcus neoformans*., *Journal of Antimicrobial Chemotherapy*. **42**: 439-444.
- (5) Linares M *et al.* (2005). Susceptibility of filamentous fungi to voriconazole tested by two microdilution methods. *Journal of Clinical Microbiology*. **43**: 250-253.
- (6) Kartsonis N *et al.* (2005). Caspofungin susceptibility testing of isolates from patients with esophageal candidiasis or invasive candidiasis: relationship of MIC to treatment outcome. *Antimicrobial Agents and Chemotherapy*. **49**: 3616-3623.
- (7) Pfaller M A *et al.* (2004). Clinical evaluation of a dried commercially prepared microdilution panel for antifungal susceptibility testing of five antifungal agents against *Candida* spp. and *Cryptococcus neoformans*. *Diagnostic Microbiology and Infectious Diseases*. **50** :113-117.
- (8) Pfaller M A *et al.* (2004). Clinical evaluation of the Sensititre YeastOne® antifungal susceptibility testing of the new triazoles voriconazole, posiconazole and ravuconazole. *Journal of Clinical Microbiology*. **42**: 4577-4580.
- (9) Espinel-Ingroff A *et al.* (2004). Multicenter comparison of the Sensititre YeastOne® colorimetric antifungal panel with the NCCLS, M27-A2 reference method for testing new antifungal agents against clinical isolates of *Candida* spp. *Journal of Clinical Microbiology* **42**: 718- 721.
- (10) Linares M *et al.* (2004). Comparison of two microdilution methods for testing susceptibility of *Candida* spp. to voriconazole. *Journal of Clinical Microbiology* **42**: 899-902.
- (11) Espinel-Ingroff A *et al.* (1999). Multicenter comparison of the Sensititre YeastOne® colorimetric antifungal panel with the National Committee for Clinical Laboratory Standards M27-A reference method for testing clinical isolates of common and emerging *Candida* spp., *Cryptococcus* spp., and other Yeasts and Yeast-Like Organisms. *Journal of Clinical Microbiology* **37**: 591-595.
- (12) Barry A L *et al.* (2000). Quality control limits for broth microdilution susceptibility tests of ten antifungal agents. *Journal of Clinical Microbiology* **38**: 3457-3459.
- (13) Canton E *et al.* (2005). Sensititre YeastOne® Caspofungin susceptibility testing of *Candida* clinical isolates: Correlation with results of NCCLS, M27-A2 multicenter study . *Antimicrobial Agents and Chemotherapy*. **49**: 1604 - 1607.
- (14) Holliday N M *et al.* (2007) A reproducibility study with two echinocandins using the Sensititre YeastOne® susceptibility plate. *American Microbiology Abstract* **C-191**.

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- (16) Espinel-Infgraff A (2006). Comparison of three commercial assays and a modified disk diffusion assay with two broth microdilution reference assays for testing Zygomycetes, *Aspergillus* spp, *Candida* spp, and *Cryptococcus neoformans* with posaconazole and Amphotericin B. *Journal of Clinical Microbiology* **44**: 3616-3622
- (17) Alexander B D *et al.* (2007). Comparative evaluation of Etest and Sensititre YeastOne® panels against the Clinical Laboratory Standards Institute M27-A2 reference broth microdilution method for testing *Candida* susceptibility to seven antifungal agents. *Journal of Clinical Microbiology* **45**: 698-706.
- (18) Eagle H. (1948). A paradoxical zone phenomenon in the bactericidal action of penicillin *in vitro*. *Science*, 1948; 107: 44-45
- (19) CLSI M27M44S. Performance Standards For Antifungal Susceptibility Testing Of Yeasts – 3rd Edition. CLSI (Clinical & Laboratory Standards Institute), 940 West Valley Road, Suite 1400, Wayne, PA 19087, US.